### The following products (code Nos.) are used in our application examples:

٠

٠

Pre-po										
0	Circular brush with Chungking bristles. Particularly beneficial are the Labor and Prima brushes with two rows of bristles,									
	converging type. (Pl	ease choo	se your fa	vorite brush from o	our catalogu	e, pages 3	7, 39 and4	1.)		
0	Circular brush – go					32/90		•	12 100)	
0	Lathe brush – goat					152/50			12 100)	
	Lathe brush – goat					152/60		(SBA 1	12 200)	
	Lathe brush – goat					162/48		•	12 100)	
0	Miniature brushes,					e brush, g				
	Ø 19 mm		(MBH 12 2		Ø 19 mm		122/19	•	12 200)	
	Ø 21 mm		(MBH 12 3		Ø 21 mm		122/21	•	12 300)	
	Ø 23 mm		(MBH 12 4	,	Ø 23 mm		122/23	•	12 400)	
	Ø 25 mm		(MBH 12 5	,	Ø 25 mm		122/25	•	12 500)	
0	Buffing wheel for pre-polishing 31/100KK Ø 100 mm B= 7.5 mm (SWK 24 100)									
	31/100LK Ø 100 mm B= 7.5 mm (SWL 24 100) Buffing wheels, cotton cloth: 29/100KK fine (SWK 22 100) 29a/100KK fine (SWK 22 200)							22 200)		
0	builling wheels, co		29/100Kr	```	,		KK coarse	•	,	
0	Lathe brushes with	n discs of		``	,	156/50			22 100)	
0	Lattie brashes with	1 41303 01		Ø 90 m		156/90		•	22 200)	
0	Miniature brushes	with discs	sof	6 50 m		e brushes	with disc	•	22 200)	
0	cotton cloth, mour				cotton cloth, unmounted					
	Ø 19 mm	102N/19	(MBH 22 2	200)	Ø 19 mm	•	129N/19	(MBU :	22 200)	
	Ø 21 mm	102N/21	•	,	Ø 21 mm		129N/21		,	
	Ø 23 mm	102N/23	•	,	Ø 23 mm		129N/23	•	,	
	Ø 25 mm	102N/25	(MBH 22 5	500)́	Ø 25 mm		129N/25	(MBU :	22 500)	
High-gl	oss polishing									
0 0	Buffing wheel, cot	ton cloth s	ee above							
0	Pleated buffing wh			29d/100KK		(SWK 28	100)			
0	Lathe brushes with	n discs of	cotton clo	oth see above		·	,			
0	Circular brushes w	ith cotton	yarn							
	20/80 KK	Ø 80 mm	-	W = approx. 16 m	m	(SRB 30	100)			
	20/90 KK	Ø 90 mm		W = approx. 16 m	m	(SRB 30 2	200)			
0	Circular brushes w		-							
	20L/80 KK	Ø 80 mm		W = approx. 13 m		(SRB 31	,			
	20L/90 KK	Ø 90 mm		W = approx. 13 m	m	(SRB 31 2	200)			
0	Circular brushes w		•	14/			000			
	20G/90 KK	Ø 90 mm		W = approx. 15 m	m	(SRB 32 2	200)			
0	Circular brushes, v 20W/80 KK	Ø 80 mm		W = approx. 16 m	m	(SRB 33	100)			
		Ø 90 mm					,			
0	20W/90 KK Ø 90 mm W = approx. 16 mm (SRB 33 200) Lathe brushes with cotton or flax yarn									
0	Cotton white			W = approx. 9 mm	'n	(SBA 30 <sup>-</sup>	100)			
	Flax yarn grey			W = approx. 9 mm		(SBA 31	,			
0	Miniature brushes,			••		<b>\</b> -	/			
0	Miniature brushes									
0	Miniature brushes,	cotton ya	rn, moun	ted	Miniature	e brushes	, cotton ya	arn, un	mounted	
	Ø 19 mm	103/19	(MBH 30 2	200)	Ø 19 mm		125/19	(MBU 3	30 200)	
	Ø 21 mm	103/21	(MBH 30 3	300)	Ø 21 mm		125/21	(MBU :	30 300)	
	Ø 23 mm		(MBH 30 4		Ø 23 mm		125/23		30 400)	
	Ø 25 mm		(MBH 30 5		Ø 25 mm		125/25		30 500)	
0	Miniature brushes,	•		•		e brushes	, grey lax 103L/19		21 200)	
	Ø 19 mm Ø 21 mm	103d/19 103d/21	(MBH 30 2	,	Ø 19 mm Ø 21 mm		103L/19 103L/21	•	31 300)	
	Ø 23 mm		(MBH 30 4		Ø 23 mm		103L/21 103L/23	•	31 400)	
	Ø 25 mm	103d/25			Ø 25 mm		103L/25	•	31 500)	
0	Miniature brushes				0 20 1111		1000/20		01 000)	
0	Miniature brushes				MB. leath	ner discs,	unmounte	ed		
-	Ø 19 mm	100/19		(MBH 00 200)	Ø 19 mm		127/19		00 200)	
	Ø 21 mm	100/21		(MBH 00 300)	Ø 21 mm		127/21	•	00 300)	
	Ø 23 mm	100/23		(MBH 00 400)	Ø 23 mm		127/23	•	00 400)	
	Ø 25 mm	100/25		(MBH 00 500)	Ø 25 mm		127/25	(MBU	00 500)	
0	Miniature brushes,		•		er), mounte	ed				
	Ø 19 mm	102NL/19		(MBH 02 200)						
	Ø 21 mm	102NL/21		(MBH 02 300)						
	Ø 23 mm	102NL/23		(MBH 02 400)						

## Wishing you good polishing results



Josef-Schüttler-Str. 49 78224 Singen, Germany

Tel: 07731-947 220 Fax: 07731-947 221



# **Polishing Acrylic Veneers**

We recommend the following articles in our range of products. In this information brochure we mainly focus on the particularly tricky places 12 - 1 - 2

for polishing.		
Circular brus Chungking bristles	shes with goat hair	Lathe brush white goat h
AND OF THE PARTY O		
But	ffing wheel ton cloth	0
	Lathe brush cotton cloth	
		Circular brushes: glazed yarn
		wool cotton yarn grey flax yarn
		Caller
Miniature brushes:		Cotton yarn
	Cotton o	cloth



nite goat hair



Miniature brush goat hair



Buffing wheel for pre-polishing





Pleated buffing wheel

> Buffing wheel combined leather/cloth







Lathe brushes with cotton yarn (white) flax yarn (grey)

Blue polishing paste





Flax yarn, grey



Chamois leather and cloth



Chamois leather

### **Pre-polishing**

Our Circular Chungking Brushes are excellently suited for pre-polishing acrylic veneers (if the surfaces are not yet smooth).

In this case we recommend using our 2-3 row Circular Labor Brushes (converging) with a low contact pressure using pumice and applying a speed of 3,000 rpm.

Example : 3/28/55 Labor, converging type (RLH 32 100)

If you have already obtained a relatively smooth surface, we recommend applying a soft Circular Goat Hair Brush at a maximum speed of 4,000 rpm. The Lathe Goat Hair Brush (maximum speed of 4,000 rpm.) is also particularly wellsuited and even enables interdental pre-polishing of the acrylic veneer. In addition, this interdental area can also be pre-polished very well manually using the Miniature Goat Hair Brush at a maximum speed of 12,000 rpm.

Example : 112/21 mm (MBH 12 300)





For this type of application we can also recommend using our special Buffing Wheel for prepolishing, a Buffing Wheel or a Lathe Brush with cotton cloth together with wet pumice.

Note: After using pumice with the buffing wheels, please make sure that they are rinsed thoroughly and that you do not exceed a speed of 4,000 rpm.

#### Do you wish for better control?

If so, we recommend using our Miniature Brush (cotton cloth) together with pumice at a maximum speed of 12,000 rom. Please make sure in this case, too, that the brush is cleaned thoroughly after use.



Example: 102N/21 mm (MBH 22 300)

### **High-gloss polishing**

In our range of products you will find a large number of brushes that are excellently suited for the high-gloss polishing of veneering resins. We are only able to give you some recommendations, information and advice for this.

You can either use the **Buffing Wheel (cotton cloth)**, the Pleated Buffing Wheel or the Lathe Brush (cotton cloth) at a maximum speed of 4,000 rpm. By applying our blue polishing paste you will obtain a beautiful gloss on your acrylic veneer.



Our circular brushes are suitable for high-gloss polishing using blue polishing paste and are made of

cotton yarn grey flax yarn glazed varn and wool

These circular brushes with a diameter of 80 or 90 mm are used particularly for polishing large acrylic veneers at a speed of 4,000 rpm.

Example : Flax yarn 20L/90 mm (SRB 31 200)

Particularly beneficial for splinted crowns and bridges with acrylic veneers, however, are the lathe brushes with cotton yarn (white)

### flax yarn (grey),

as they are well-suited for polishing the interdental spaces. We recommend using blue polishing paste at a speed of 3,000 rpm.

Example : Cotton yarn 155/50 mm (SBA 30 100)

#### Now and again, there are still a few small places that have not been optimally polished – are you familiar with this problem?

You will find the solution with our Miniature Brushes (MB).

If there is still an interdental place that still needs some more polishing, we again recommend using our Goat Hair White. In contrast to pre-polishing, apply the **blue polishing paste (P3)** here. In this case, the maximum speed of 12,000 rpm should be reduced, if necessary, as the removal rate of this brush is often underestimated. Please also make sure to use a lower contact pressure and to choose a brush with a larger diameter (21 mm or more) to prevent any risk of overheating.

For this application our miniature brushes with **cotton** or **flax yarn**, or cotton cloth or flannel are also ideal, and are used with the hand-piece at a maximum speed of 12,000 rpm.

#### Have you ever tried leather?

Our miniature brushes with chamois leather, or a combination of leather/cloth/leather produce a particularly high gloss at a maximum speed of 12,000 rpm. When using the Miniature Brush (leather) you do not need any additional polishing paste. However, vou should be fully aware of the risk of overheating. The combined leather/cloth/leather variant generates lower heat than the leather brush without cloth. In addition, you can apply **blue** polishing paste to the cotton cloth.

Example : Buffing wheel, cotton cloth, 29/100 mm (SWK 22 100)







Goat hair 112/21 mm (MBH 12 300)



Cotton 103/21 mm (MBH 30 300)



Leather/cloth 102NL/21 mm (MBH 02 300)